

# **Long-read amplicon sequencing for microbiome analysis**

**Yan Hui**

**Postdoc**

**Department of Food Science**

**E-mail: [huiyan@food.ku.dk](mailto:huiyan@food.ku.dk)**

UNIVERSITY OF COPENHAGEN



# Table of Contents

- **Reference-based and reference-free strategies for lengthy amplicon analysis**
- **De novo OTU picking from long amplicons with LACA**
- **Use NART for long amplicon profiling by read classification**
- **Exercises**



# Reference-dependent vs Reference-free analysis

OUTPUT	RRF-dependent	REF-free
Representative sequences	No	Yes
Phylogenetic tree	No	Yes



## EPI2ME

Per-read query against a known database:

- Limited by database
- No OTUs

## illumina<sup>®</sup>

The consensus from built clusters:


- Clustering by identity, etc.

# Sequencing errors

Phred quality scores  $Q$  are logarithmically related to the base-calling error probabilities  $P$  and defined as

$$Q = -10 \log_{10} P.$$


illumina®




**PCR >>>  
More errors in the  
end**

**High error rate**

R9.4.1



R10.3




**Bio-pore >>>  
Random error**

# Molecule-level correction

## ONT Duplex

Linear dsDNA molecule adapted on both ends and first strand sequenced

Second strand captured and sequenced subsequently



# PacBio Circular Consensus Sequencing (CCS)

Start with high-quality double stranded DNA


Ligate SMRTbell adapters and size select

Anneal primers and bind DNA polymerase

Circularized DNA is sequenced in repeated passes


The polymerase reads are trimmed of adapters to yield subreads

Consensus is called from subreads




# Molecule-level correction

## Specialized library preparation with e.g., UMI



# Clustering-based correction



Metagenomics: microbes in uneven abundance

UMI -> different template (including the phylogenetically same one)

# Clustering-based correction

## Troubles with long-read alignment

- **Pairwise alignment**

Time complexity

$$O(C(n,2)/2) = O(n!(n-1)!/2)$$


For amplicons,  $n$  can be millions if reads are pooled.

- **Noisy alignment with long reads**

The relatively **high** error rate in relatively **long** reads




# K-mers binning





	atcga	tcgac	...	cgaaa
read1	1	6		7
read2	3	1		3
...				
readn	5	2		2

- Computers prefer k-mers than text: blast, binning
- Unique k-mer patterns between genomes



# Pre-cluster: Use 5-kmer profiles to bin ONT reads



```

AAAGCCATATTATTCTTATGCTCTGGGTCTATTATTACACA-GC
AAAG-ATATTATTTCTTTTGTTCAGGTCAATATTCTTCATGC
AAAGCAATGCTGTTCTATGCTCCGGATCAATTATTCACA-GT
AAAGCCATACTATTTATATGCTCCGGGTCCATCATCCACA-AC
AAAGCCATACTATTCATATGCTCAGGGTCAATTATCCACA-GC
AAAGCTATACTATTCCTATGCTCAGGGTCCATTATTCATA-GC
AAAGCTATACTCTTCTTATGCTCTGGATCCATCATCCACA-AC
AAAGCCATATTATTCCTATGCTCGGGTGCATTATCCACA-AC
AAAGCCATACTATTCCTATGCTCCGGTTCATTATCCACA-AT

```

	atcga	tcgac	...	cgaaa
read1	1	10		3
read2	1	10		4
...				
readn	1	10		3

Take **blast** result as an example

Raw reads within the cluster

Raw read

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
✓	<a href="#">Limosilactobacillus fermentum strain 9-4 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1827	9070	99%	0.0	90.92%	2085632	<a href="#">CP076082.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HFD1 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9096	99%	0.0	90.86%	2101878	<a href="#">CP050919.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain AGR1487 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9113	99%	0.0	90.86%	1939032	<a href="#">CP047585.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain USM 8633 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9085	99%	0.0	90.86%	2238401	<a href="#">CP045034.1</a>
✓	<a href="#">Lactobacillus fermentum strain SL1-1 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1513	<a href="#">MN435796.1</a>
✓	<a href="#">Lactobacillus fermentum strain IITKGP-BT13 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1513	<a href="#">MN267492.1</a>
✓	<a href="#">Lactobacillus fermentum strain BioE LF11 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1512	<a href="#">MK779053.1</a>
	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
✓	<a href="#">Limosilactobacillus fermentum strain B1 28 chromosome</a>	<a href="#">Limosilacto...</a>	1884	9126	100%	0.0	89.82%	1905587	<a href="#">CP039750.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HBUAS54312 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	89.22%	1498	<a href="#">MH817761.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HBUAS62516 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	89.22%	1498	<a href="#">ON005289.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HFD1 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9039	100%	0.0	89.16%	2101878	<a href="#">CP050919.1</a>
✓	<a href="#">Limosilactobacillus fermentum 3872 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9033	100%	0.0	89.16%	2297851	<a href="#">CP011536.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain ACA-DC 179 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9022	100%	0.0	89.16%	2149913	<a href="#">CP082359.1</a>

Take **blast** result as an example

## Denoised consensus

## Polished consensus

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
✓	<a href="#">Limosilactobacillus fermentum strain AGR1485 chromosome...</a>	<a href="#">Limosila...</a>	2728	13615	100%	0.0	100.00%	2226862	<a href="#">CP047584.1</a>
✓	<a href="#">Lactobacillus fermentum strain shebah-101 16S ribosomal R...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1494	<a href="#">MN625236.1</a>
✓	<a href="#">Lactobacillus fermentum strain HB 16S ribosomal RNA gene...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1509	<a href="#">MN589591.1</a>
✓	<a href="#">Lactobacillus fermentum strain SL5-1 16S ribosomal RNA ge...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1513	<a href="#">MN435802.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain B1 28 chromosome</a>	<a href="#">Limosila...</a>	2728	13574	100%	0.0	100.00%	1905587	<a href="#">CP039750.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HDB1096 16S ribosom...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1492	<a href="#">MK537375.1</a>
✓	<a href="#">Lactobacillus fermentum strain LF 16S ribosomal RNA gene...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1564	<a href="#">MK245999.1</a>
✓	<a href="#">Lactobacillus fermentum strain LMEM36 16S ribosomal RNA...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1545	<a href="#">MK239985.1</a>
✓	<a href="#">Lactobacillus fermentum strain LMEM19 16S ribosomal RNA...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1529	<a href="#">MK239955.1</a>
✓	<a href="#">Lactobacillus fermentum strain S1 16S ribosomal RNA gene...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1531	<a href="#">MK226442.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain MTCC 5898 chromosome</a>	<a href="#">Limosila...</a>	2728	13600	100%	0.0	100.00%	2098685	<a href="#">CP035904.1</a>
✓	<a href="#">Lactobacillus fermentum strain LMEM 5 16S ribosomal RNA...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1528	<a href="#">MK418591.1</a>
✓	<a href="#">Lactobacillus fermentum strain LMEM 37 16S ribosomal RN...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1557	<a href="#">MK418588.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain LDTM 7301 chromoso...</a>	<a href="#">Limosila...</a>	2728	13593	100%	0.0	100.00%	2046196	<a href="#">CP031195.1</a>
✓	<a href="#">Lactobacillus fermentum strain PRS1 16S ribosomal RNA ge...</a>	<a href="#">Limosila...</a>	2728	2728	100%	0.0	100.00%	1515	<a href="#">MH472943.1</a>

# *De novo* OTU picking from long amplicons with **LACA**

UNIVERSITY OF COPENHAGEN



# LACA: an automatic workflow for Long Amplicon Consensus Analysis



- GitHub: <https://github.com/yanhui09/laca>

## Example

```
laca init -b /path/to/basecalled_fastqs -d /path/to/database    # init config file and check
laca run all                                           # start analysis
```





# Use **NART** for long amplicon profiling by read classification

UNIVERSITY OF COPENHAGEN






# NART: A tool for Nanopore Amplicon Real-Time analysis



- GitHub: <https://github.com/yanhui09/nart>
- Demo video:  
<https://www.youtube.com/watch?v=TkdJGLOscPg>



# Directed Acyclic Graph (DAG)

**a. emu****b. minimap2lca****c. blast2lca**

# Lowest Common Ancestor by read classification (minimap2lca, blast2lca)


*Limosilactobacillus fermentum*

Raw read


	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
✓	<a href="#">Limosilactobacillus fermentum strain 9-4 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1827	9070	99%	0.0	90.92%	2085632	<a href="#">CP076082.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HFD1 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9096	99%	0.0	90.86%	2101878	<a href="#">CP050919.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain AGR1487 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9113	99%	0.0	90.86%	1939032	<a href="#">CP047585.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain USM 8633 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1823	9085	99%	0.0	90.86%	2238401	<a href="#">CP045034.1</a>
✓	<a href="#">Lactobacillus fermentum strain SL1-1 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1513	<a href="#">MN435796.1</a>
✓	<a href="#">Lactobacillus fermentum strain IITKGP-BT13 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1513	<a href="#">MN267492.1</a>
✓	<a href="#">Lactobacillus fermentum strain BioE LF11 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	90.86%	1512	<a href="#">MK779053.1</a>
	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
✓	<a href="#">Limosilactobacillus fermentum strain B1 28 chromosome</a>	<a href="#">Limosilacto...</a>	1884	9126	100%	0.0	89.82%	1905587	<a href="#">CP039750.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HBUAS54312 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	89.22%	1498	<a href="#">MH817761.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HBUAS62516 16S ribosomal RNA gene, partial sequence</a>	<a href="#">Limosilacto...</a>	1823	1823	99%	0.0	89.22%	1498	<a href="#">ON005289.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain HFD1 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9039	100%	0.0	89.16%	2101878	<a href="#">CP050919.1</a>
✓	<a href="#">Limosilactobacillus fermentum 3872 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9033	100%	0.0	89.16%	2297851	<a href="#">CP011536.1</a>
✓	<a href="#">Limosilactobacillus fermentum strain ACA-DC 179 chromosome, complete genome</a>	<a href="#">Limosilacto...</a>	1820	9022	100%	0.0	89.16%	2149913	<a href="#">CP082359.1</a>

# Emu: Species-level abundance estimation through an expectation–maximization algorithm


## a Maximum likelihood




## b Expectation maximization




Map reads ( $R$ ) to database sequences ( $S$ )



Alignment type probabilities [ $P(c)$ ]




Alignment probabilities [ $P(r/t)$ ]



Expectation

Species-level probabilities [ $P(t/t)$ ]




Composition estimate [ $F(t)$ ]



Maximization

Trim noise and output the result



# nart & nawf



**NART** is composed of two sets of scripts: **nart** and **nawf**, which controls real-time analysis and workflow performance, respectively.

Usage: nart [OPTIONS] COMMAND [ARGS]...

NART: A tool for Nanopore Amplicon Real-Time (NART) analysis. To follow updates and report issues, see: <https://github.com/yanhui09/nart>.

## Options:

-v, --version Show the version and exit.  
-h, --help Show this message and exit.

## Commands:

monitor Start NART to monitor a directory.  
run Start NART workflow.  
visual Start NART app to interactively visualize the results.

Usage: nawf [OPTIONS] COMMAND [ARGS]...

NAWF: A sub-tool to run Nanopore Amplicon WorkFlow. The workflow command initiates the NAWF in a single batch, using either a fastq file from one ONT run or a fastq file generated during sequencing. To follow updates and report issues, see: <https://github.com/yanhui09/nart>.

## Options:

-v, --version Show the version and exit.  
-h, --help Show this message and exit.

## Commands:

config Generate the workflow config file.  
run Start workflow in a single batch.

# Usage



## Amplicon analysis in single batch

`nawf` can be used to profile any single basecalled `fastq` file from a Nanopore run or batch.

```
nawf config -b /path/to/basecall_fastq -d /path/to/database # init config file and check
nawf run all # start analysis
```



## Real-time analysis

`nart` provide utils to record, process and profile the continuously generated `fastq` batch.

Before starting real-time analysis, you need `nawf` to configure the workflow according to your needs.

```
nawf config -d /path/to/database # init config file and check
```



In common cases, you need three independent sessions to handle monitor, process and visualization, respectively.

1. Monitor the basecall output and record

```
nart monitor -q /path/to/basecall_fastq_dir # monitor basecall output
```



2. Start amplicon analysis for new fastq

```
nart run -t 10 # real-time process in batches
```



3. Update the feature table for interactively visualize in the browser

```
nart visual # interactive visualization
```



# RT-philosophy



## ONT sequencing and basecalling in batches

- nart monitor => fqs.txt (record fastq files)
- nart run => nawf (start the workflow in batches & update the feature table)
- nart visual => interactively visualize profiles.

# Exercises

UNIVERSITY OF COPENHAGEN



# Exercise

- MAC2023:

<https://yanhui09.github.io/MAC2023/>



**Cross-platform  
support, incl. MacOS**

**CONDA®**



**Linux/amd64 platform**



# Thanks